

# TQS

Tétanos Quick Stick

SINGLE TEST

Save time, 1 test for 1 patient !



Get the immune status in 10 min<sup>3</sup>

## UNDERIMMUNIZATION

Increased risk to contract tetanus



## OVERIMMUNIZATION

Increased costs and side effects

	Whole blood	Plasma	Serum
<b>Sensitivity (VPP)</b>	97.0 % n=169 (100 %)	97.2 % n=216 (100 %)	98.6 % n=213 (100 %)
<b>Specificity (VPN)</b>	100 % n=75 (93.8 %)	100 % n=43 (87.8 %)	100 % n=44 (93.6 %)

1. D.T., C.P. 2009. Rapid immunochromatographic tests. High Authority of Health setting evidence of antitetanic immunoprotection in emergency settings. P36. [19/04/2019] <https://docplayer.fr/863267-Mise-en-evidence-de-l-immunoprotection-antitetanique-en-contexte-d-urgence-evaluation-des-tests-rapides-immunochromatographiques.html>
2. N'Diaye et al. BMC Infectious Diseases 2014, 14:603 <http://www.biomedcentral.com/1471-2334/14/603>
3. For more information, please refer to the paragraph «Result interpretation» of the technical sheet.

**Eurobio Scientific**

7 Avenue de Scandinavie ZA de Courtaboeuf

91 953 LES ULIS CEDEX B France

Phone : +33 1 69 79 64 80 Fax : +33 1 69 79 05 35

marketing@eurobio-scientific.com www.eurobio-scientific.com

EUR0719-TQS03-EN

Eurobio Scientific (France)



This medical device for *in vitro* diagnosis is a regulated health product and thus bears the CE mark.



SEE performance study as mentioned in the technical sheet.